

## Product datasheet for TA500382AM

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### SIGLEC9 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI6D8]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI6D8

**Applications:** IF, IHC, WB

Recommended Dilution: WB 1:500~2000, IHC 1:50, IF 1:50

Reactivity: Human, Dog, Monkey

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human SIGLEC9 (NP\_055256) produced in

HEK293T cell.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 50.1 kDa

**Gene Name:** sialic acid binding lg like lectin 9

Database Link: NP 055256

Entrez Gene 719831 MonkeyEntrez Gene 27180 Human

Q9Y336

**Background:** Putative adhesion molecule that mediates sialic-acid dependent binding to cells.Preferentially

binds to alpha-2,3- or alpha-2,6-linked sialic acid. The sialic acid recognition site may be

masked by cis interactions with sialic acids on the same cell surface.


Synonyms: CD329; CDw329; FOAP-9; OBBP-LIKE; siglec-9

**Protein Families:** Druggable Genome, Transmembrane







# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SIGLEC9 ([RC206674], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SIGLEC9. Positive lysates [LY402335] (100ug) and [LC402335] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-SIGLEC9 monoclonal antibody.



Immunohistochemical staining of paraffinembedded Human breast tissue within the normal limits using anti-SIGLEC9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500382])





Immunofluorescent staining of HT29 cells using anti-SIGLEC9 mouse monoclonal antibody ([TA500382]).