

### Product datasheet for TA500447BM

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## B Raf (BRAF) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1G8]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1G8

Applications: IHC, WB

Recommended Dilution: WB 1:2000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human BRAF expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 84.4 kDa

**Gene Name:** B-Raf proto-oncogene, serine/threonine kinase

Database Link: NP 004324

Entrez Gene 109880 MouseEntrez Gene 114486 RatEntrez Gene 673 Human

P15056





Background:

This gene encodes a protein belonging to the raf/mil family of serine/threonine protein kinases. This protein plays a role in regulating the MAP kinase/ERKs signaling pathway, which affects cell division, differentiation, and secretion. Mutations in this gene are associated with cardiofaciocutaneous syndrome, a disease characterized by heart defects, mental retardation and a distinctive facial appearance. Mutations in this gene have also been associated with various cancers, including non-Hodgkin lymphoma, colorectal cancer, malignant melanoma, thyroid carcinoma, non-small cell lung carcinoma, and adenocarcinoma of lung. A pseudogene, which is located on chromosome X, has been identified for this gene. [provided by RefSeq]


Synonyms: B-raf; B-RAF1; BRAF1; NS7; RAFB1

**Protein Families:** Druggable Genome, Protein Kinase

Acute myeloid leukemia, Bladder cancer, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Focal adhesion, Glioma, Insulin signaling pathway, Long-term depression, Long-term potentiation, MAPK signaling pathway, Melanoma, mTOR signaling pathway, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer, Regulation of actin cytoskeleton, Renal cell carcinoma, Thyroid cancer, Vascular smooth muscle contraction


# **Product images:**

**Protein Pathways:** 




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY BRAF ([RC211013], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BRAF (1:500).






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-BRAF mouse monoclonal antibody. ([TA500447])




Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-BRAF mouse monoclonal antibody. ([TA500447])




Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-BRAF mouse monoclonal antibody. ([TA500447])





Immunohistochemical staining of paraffinembedded Human bladder tissue within the normal limits using anti-BRAF mouse monoclonal antibody. ([TA500447])



Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-BRAF mouse monoclonal antibody. ([TA500447])