

Product datasheet for TA500494BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

HSP90AB1 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4C10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4C10

Applications: WB

Recommended Dilution: WB 1:2000

Reactivity: Human, Mouse, Rat, Dog, Monkey

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human HSP90AB1 expression

vector

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 83.3 kDa

Gene Name: heat shock protein 90 alpha family class B member 1

Database Link: NP 031381

Entrez Gene 15516 MouseEntrez Gene 301252 RatEntrez Gene 474919 DogEntrez Gene

702293 MonkeyEntrez Gene 3326 Human

P08238





Background: HSP90 proteins are highly conserved molecular chaperones that have key roles in signal

transduction, protein folding, protein degradation, and morphologic evolution. HSP90 proteins normally associate with other cochaperones and play important roles in folding newly synthesized proteins or stabilizing and refolding denatured proteins after stress. There are 2 major cytosolic HSP90 proteins, HSP90AA1 (MIM 140571), an inducible form, and HSP90AB1, a constitutive form. Other HSP90 proteins are found in endoplasmic reticulum

(HSP90B1; MIM 191175) and mitochondria (TRAP1; MIM 606219)


Synonyms: D6S182; HSP84; HSP90B; HSPC2; HSPCB

Protein Families: Druggable Genome, Stem cell - Pluripotency

Protein Pathways: Antigen processing and presentation, NOD-like receptor signaling pathway, Pathways in

cancer, Progesterone-mediated oocyte maturation, Prostate cancer

Product images:



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-HSP90AB1 monoclonal antibody.