

## **Product datasheet for TA500496S**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### STK39 Mouse Monoclonal Antibody [Clone ID: OTI4H3]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4H3

**Applications:** FC, IF, IHC, WB

**Recommended Dilution:** WB 1:1000~2000, IHC 1:50, IF 1:50, FLOW 1:100

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human STK39 (NP\_037365) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 59.3 kDa

**Gene Name:** serine/threonine kinase 39

Database Link: NP 037365

Entrez Gene 53416 MouseEntrez Gene 54348 RatEntrez Gene 27347 Human

Q9UEW8

**Background:** This gene encodes a serine/threonine kinase that is thought to function in the cellular stress

response pathway. The kinase is activated in response to hypotonic stress, leading to phosphorylation of several cation-chloride-coupled cotransporters. The catalytically active kinase specifically activates the p38 MAP kinase pathway, and its interaction with p38 decreases upon cellular stress, suggesting that this kinase may serve as an intermediate in

the response to cellular stress. [provided by RefSeq]







**Synonyms:** DCHT; PASK; SPAK

**Protein Families:** Druggable Genome, Protein Kinase


# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY STK39 ([RC223981], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-STK39. Positive lysates [LY402227] (100ug) and [LC402227] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by usin g anti-STK39 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).




Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-STK39 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500496])





Immunofluorescent staining of COS7 cells using anti-STK39 mouse monoclonal antibody ([TA500496]).



HEK293T cells transfected with either [RC223981] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-STK39 antibody ([TA500496]), and then analyzed by flow cytometry.