

Product datasheet for TA500536BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

TrkC (NTRK3) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI24B3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI24B3
Applications: IHC, WB

Reactivity: WB 1:2000, IHC 1:50 Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human NTRK3 expression vector

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 94.4 kDa

Gene Name: neurotrophic receptor tyrosine kinase 3

Database Link: NP 001007157

Entrez Gene 18213 MouseEntrez Gene 29613 RatEntrez Gene 4916 Human

Q16288

Background: This gene encodes a member of the neurotrophic tyrosine receptor kinase (NTRK) family. This

kinase is a membrane-bound receptor that, upon neurotrophin binding, phosphorylates itself and members of the MAPK pathway. Signalling through this kinase leads to cell differentiation and may play a role in the development of proprioceptive neurons that sense body position. Mutations in this gene have been associated with medulloblastomas, secretory breast

carcinomas and other cancers.






Synonyms: gp145(trkC); TRKC


Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Protein Pathways: Neurotrophin signaling pathway


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NTRK3 ([RC203781], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NTRK3.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-NTRK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500536])




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-NTRK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500536])






Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-NTRK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500536])




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-NTRK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500536])




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-NTRK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500536])






Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-NTRK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500536])




Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-NTRK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500536])



Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-NTRK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500536])





Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-NTRK3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500536])