

### Product datasheet for TA500560BM

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Protein kinase Y linked (PRKY) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4G7]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4G7

Applications: IF, IHC, WB

Recommended Dilution: WB 1:500, IHC 1:50, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full-length protein expressed in 293T cell transfected with human PRKY expression vector

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

Predicted Protein Size: 31.7 kDa

**Gene Name:** protein kinase Y-linked (pseudogene)

Database Link: NP 002751

Entrez Gene 5616 Human

<u>O43930</u>

Background: This gene is similar to the protein kinase, X-linked gene in the pseudoautosomal region of the

X chromsoome. The gene is classified as a transcribed pseudogene because it has lost a coding exon that results in all transcripts being candidates for nonsense-mediated decay (NMD) and unlikely to express a protein. Abnormal recombination between this gene and a

related gene on chromosome X is a frequent cause of XX males and XY females.







Synonyms: OTTHUMP00000033227; protein kinase, Y-linked

**Protein Families:** Druggable Genome, Protein Kinase

## **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PRKY ([RC210438], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PRKY. Positive lysates [LY400977] (100ug) and [LC400977] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of breast tissue using anti-PRKY mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)




Immunohistochemical staining of paraffinembedded colon tissue within the normal limits using anti-PRKYmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)



Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-PRKYmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)




Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-PRKYmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)




Immunohistochemical staining of paraffinembedded lung tissue within the normal limits using anti-PRKYmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)






Immunohistochemical staining of paraffinembedded Ovary tissue within the normal limits using anti-PRKYmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-PRKYmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)




Immunohistochemical staining of paraffinembedded pancreas tissue within the normal limits using anti-PRKYmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)






Immunohistochemical staining of paraffinembedded Carcinoma of thyroid tissue using anti-PRKYmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)



Immunohistochemical staining of paraffinembedded endometrium tissue within the normal limits using anti-PRKYmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)




Immunohistochemical staining of paraffinembedded prostate tissue within the normal limits using anti-PRKYmouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)





Immunohistochemical staining of paraffinembedded Carcinoma of prostate tissue using anti-PRKYmouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500560], Dilution 1:50)



Anti-PRKY mouse monoclonal antibody ([TA500560]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PRKY ([RC210438]).