

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA500620BM

SSX2 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4D10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI4D10
Applications:	IF, IHC, WB
Recommended Dilution:	WB 1:2000, IHC 1:50, IF 1:100
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Full-length protein expressed in 293T cell transfected with human SSX2 expression vector
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	21.6 kDa
Gene Name:	SSX family member 2
Database Link:	<u>NP_003138</u> <u>Entrez Gene 6757 Human</u> <u>Q16385</u>




SSX2 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4D10] – TA500620BM

Background:The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX)
breakpoint proteins. These proteins may function as transcriptional repressors. They are also
capable of eliciting spontaneously humoral and cellular immune responses in cancer
patients, and are potentially useful targets in cancer vaccine-based immunotherapy. SSX1,
SSX2 and SSX4 genes have been involved in the t(X;18) translocation characteristically found
in all synovial sarcomas. This translocation results in the fusion of the synovial sarcoma
translocation gene on chromosome 18 to one of the SSX genes on chromosome X. The
encoded hybrid proteins are probably responsible for transforming activity. Two transcript
variants encoding distinct isoforms have been identified for this gene.


Synonyms: CT5.2; CT5.2A; HD21; HOM-MEL-40; SSX

Protein Families: Druggable Genome, Transcription Factors


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY SSX2 (Cat# [RC201214], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SSX2(Cat# [TA500620]). Positive lysates [LY418870] (100ug) and [LC418870] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 8 different cell lines by using anti-SSX2 monoclonal antibody at 1:200 dilution.



Immunohistochemical staining of paraffinembedded Adenocarcinoma of breast tissue using anti-SSX2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)

Immunohistochemical staining of paraffinembedded colon tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)


Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-SSX2mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)



Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)

Immunohistochemical staining of paraffinembedded Ovary tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)


Immunohistochemical staining of paraffinembedded pancreas tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)



Immunohistochemical staining of paraffinembedded Carcinoma of thyroid tissue using anti-SSX2mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)

Immunohistochemical staining of paraffinembedded endometrium tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)

Immunohistochemical staining of paraffinembedded prostate tissue within the normal limits using anti-SSX2mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)



Immunohistochemical staining of paraffinembedded Carcinoma of bladder tissue using anti-SSX2mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500620], Dilution 1:50)

Anti-SSX2 mouse monoclonal antibody ([TA500620]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY SSX2 ([RC201214]).