

Product datasheet for TA500685

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

BCKDH kinase (BCKDK) Mouse Monoclonal Antibody [Clone ID: OTI11C9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI11C9
Applications: IF, IHC, WB

Recommended Dilution: WB 1:2000, IHC 1:50, IF 1:100

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human BCKDK (NP_005872) produced in HEK293T

cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.8 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 46.4 kDa

Gene Name: branched chain keto acid dehydrogenase kinase

Database Link: NP 005872

Entrez Gene 12041 MouseEntrez Gene 29603 RatEntrez Gene 10295 Human

<u>O14874</u>

Background: Catalyzes the phosphorylation and inactivation of the branched-chain alpha-ketoacid

dehydrogenase complex,the key regulatory enzyme of the valine, leucine and isoleucine catabolic pathways. Key enzyme that regulate the activity state of the BCKD complex .


Synonyms: BCKDKD; BDK

Protein Families: Druggable Genome, Protein Kinase







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY BCKDK (Cat# [RC203601], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BCKDK(Cat# TA500685). Positive lysates [LY401780] (100ug) and [LC401780] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human Ovary tissue within the normal limits using anti-BCKDK mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500685)



Immunohistochemical staining of paraffinembedded Human endometrium tissue within the normal limits using anti-BCKDK mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA500685)





Anti-BCKDK mouse monoclonal antibody (TA500685) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY BCKDK ([RC203601]).