

# **Product datasheet for TA500918AM**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## HIPPI (IFT57) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI4F6]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4F6

**Applications:** IF, IHC, WB

Recommended Dilution: WB 1:2000, IHC 1:50, IF 1:100

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human IFT57(NP\_060480) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 49.1 kDa

**Gene Name:** intraflagellar transport 57

Database Link: NP 060480

Entrez Gene 73916 MouseEntrez Gene 303968 RatEntrez Gene 55081 Human

O9NWB7

**Background:** Required for the formation of cilia. Plays an indirect role in sonic hedgehog signaling, cilia

being required for all activity of the hedgehog pathway (By similarity). Has pro-apoptotic function via its interaction with HIP1,leading to recruit caspase-8 (CASP8) and trigger apoptosis. Has the ability to bind DNA sequence motif 5'-AAAGACATG-3'present in the

promoter of caspase genes such as CASP1, CASP8 and CASP10, suggesting that it may act as a

transcription regulator; however the relevance of such function remains unclear







**Synonyms:** ESRRBL1; HIPPI; MHS4R2

Protein Families: Druggable Genome
Protein Pathways: Huntington's disease

## **Product images:**







HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY IFT57 (Cat# [RC204116], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-IFT57 (Cat# [TA500918]). Positive lysates [LY402637] (100ug) and [LC402637] (20ug) can be purchased separately from OriGene.

Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-IFT57 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human) (1:200).






Immunohistochemical staining of paraffinembedded breast tissue within the normal limits using anti-IFT57 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500918], Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of colon tissue using anti-IFT57mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500918], Dilution 1:50)




Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-IFT57mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500918], Dilution 1:50)





Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-IFT57mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500918], Dilution 1:50)



Anti-IFT57 mouse monoclonal antibody ([TA500918]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY IFT57 ([RC204116]).