

## Product datasheet for TA500956AM

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Natriuretic Peptide Receptor C (NPR3) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI2G2]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI2G2

**Applications:** FC, IF, IHC, WB

Recommended Dilution: WB 1:2000, IHC 1:50, IF 1:100

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NPR3 (NP\_000899) produced in HEK293T

cell.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

Storage: Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 59.8 kDa

**Gene Name:** natriuretic peptide receptor 3

Database Link: NP 000899

Entrez Gene 18162 MouseEntrez Gene 25339 RatEntrez Gene 4883 Human

P17342

Background: The family of natriuretic peptides (see MIM 108780) elicit a number of vascular, renal, and

endocrine effects that are important in the maintenance of blood pressure and extracellular fluid volume. These effects are mediated by specific binding of the peptides to cell surface

receptors in the vasculature, kidney, adrenal, and brain






# Natriuretic Peptide Receptor C (NPR3) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI2G2] – TA500956AM

Synonyms: ANP-C; ANPR-C; ANPRC; C5orf23; GUCY2B; NPR-C; NPRC

**Protein Families:** Druggable Genome, Transmembrane

### **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NPR3 ([RC219453], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NPR3. Positive lysates [LY424462] (100ug) and [LC424462] (20ug) can be purchased separately from OriGene.






Figure from citation: Western Blot of NPR3 protein level by using anti-NPR3 antibody in HEK293 cells. Dilution: 1:1000 <u>View Citation</u>




Immunohistochemical staining of paraffinembedded Kidney tissue within the normal limits using anti-NPR3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500956], Dilution 1:50)






Immunohistochemical staining of paraffinembedded liver tissue within the normal limits using anti-NPR3mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500956], Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of ovary tissue using anti-NPR3mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA500956], Dilution 1:50)



Anti-NPR3 mouse monoclonal antibody ([TA500956]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY NPR3 ([RC219453]).





HEK293T cells transfected with either pCMV6-ENTRY NPR3 ([RC219453]) (Red) or empty vector control plasmid (Blue) were immunostained with anti-NPR3 mouse monoclonal ([TA500956]), and then analyzed by flow cytometry.