

Product datasheet for **TA500975S**


FAHD2A Mouse Monoclonal Antibody [Clone ID: OTI9G1]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI9G1
Applications:	FC, IF
Recommended Dilution:	IF 1:100, FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FAHD2A (NP_057128) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.8 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	34.6 kDa
Gene Name:	fumarylacetoacetate hydrolase domain containing 2A
Database Link:	NP_057128 Entrez Gene 51011 Human Q96GK7
Background:	May have hydrolase activity
Synonyms:	CGI-105



[View online »](#)

Product images:

Anti-FAHD2A mouse monoclonal antibody ([TA500975]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FAHD2A ([RC211128]).



HEK293T cells transfected with either [RC211128] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-FAHD2A antibody ([TA500975]), and then analyzed by flow cytometry.