

## Product datasheet for TA501164BM

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CD147 (BSG) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI9H5]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI9H5

**Applications:** FC, IF, IHC, WB

**Recommended Dilution:** WB 1:1500~3000, IHC 1:50, IF 1:100, FLOW 1:100

Reactivity: Human, Dog, Rat, Monkey

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human BSG(NP\_001719) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 40.2 kDa

**Gene Name:** basigin (Ok blood group)

Database Link: NP 001719

Entrez Gene 25246 RatEntrez Gene 476758 DogEntrez Gene 721068 MonkeyEntrez Gene 682

<u>Human</u> P35613

**Background:** The protein encoded by this gene is a plasma membrane protein that is important in

spermatogenesis, embryo implantation, neural network formation, and tumor progression.

The encoded protein is also a member of the immunoglobulin superfamily. Multiple transcript variants encoding different isoforms have been found for this gene.







Synonyms: 5F7; CD147; EMMPRIN; EMPRIN; HAb18G; OK; SLC7A11; TCSF

**Protein Families:** Druggable Genome, ES Cell Differentiation/IPS, Transmembrane


## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY BSG ([RC219464], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BSG. Positive lysates [LY419777] (100ug) and [LC419777] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-BSG monoclonal antibody.




Western blot analysis of extracts (10ug) from 10 Human tissue by using anti-BSG monoclonal antibody at 1:500 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: colon;10: spleen).






Equivalent amounts of cell lysates (10 ug per lane) of wild-type 293T cells (WT, Cat# LC810293T) and BSG-Knockout 293T cells (KO, Cat#[LC810433]) were separated by SDS-PAGE and immunoblotted with anti-BSG monoclonal antibody [TA501164] (1:2000`). Then the blotted membrane was stripped and reprobed with anti-PCNA antibody as a loading control.




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-BSG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA501164], Dilution 1:50)




Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-BSG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA501164], Dilution 1:50)






Anti-BSG mouse monoclonal antibody ([TA501164]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY BSG ([RC219464]).



Immunofluorescent staining of HT29 cells using anti-BSG mouse monoclonal antibody ([TA501164]).



HEK293T cells transfected with either [RC219464] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-BSG antibody ([TA501164]), and then analyzed by flow cytometry.