

## Product datasheet for **TA501193AM**

### IGBP1 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI4F6]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI4F6
Applications:	FC, IF, WB
Recommended Dilution:	IF 1:100, FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human IGBP1 (NP_001542) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	39.0 kDa
Gene Name:	immunoglobulin binding protein 1
Database Link:	<a href="#">NP_001542</a> <a href="#">Entrez Gene 3476 Human</a> <a href="#">P78318</a>
Background:	The proliferation and differentiation of B cells is dependent upon a B-cell antigen receptor (BCR) complex. Binding of antigens to specific B-cell receptors results in a tyrosine phosphorylation reaction through the BCR complex and leads to multiple signal transduction pathways.
Synonyms:	ALPHA-4; IBP1




[View online »](#)

Protein Families: Druggable Genome


**Product images:**



Western blot analysis of extracts (10ug) from 5 different cell lines by using anti-IGBP1 monoclonal antibody (1:200).



Anti-IGBP1 mouse monoclonal antibody ([TA501193]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY IGBP1 ([RC200387]).



HEK293T cells transfected with either [RC200387] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-IGBP1 antibody ([TA501193]), and then analyzed by flow cytometry.