

## Product datasheet for TA501206BM

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### PDE4A Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI5C6]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI5C6

**Applications:** FC, IF, IHC, WB

**Recommended Dilution:** WB 1:2000, IHC 1:50, IF 1:100, FLOW 1:100

Reactivity: Human, Rat

Host: Mouse lgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human PDE4A (NP\_006193) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 72.0 kDa

**Gene Name:** phosphodiesterase 4A

Database Link: NP 006193

Entrez Gene 5141 Human

P27815

**Background:** Cyclic nucleotides are important second messengers that regulate and mediate a number of

cellular responses to extracellular signals, such as hormones, light, and neurotransmitters. Cyclic nucleotide phosphodiesterases (PDEs) regulate the cellular concentrations of cyclic nucleotides and thereby play a role in signal transduction. PDE4A is a class IV cAMP-specific

PDE






Synonyms: DPDE2; PDE4; PDE46


Protein Families: Druggable Genome

**Protein Pathways:** Progesterone-mediated oocyte maturation, Purine metabolism


# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PDE4A ([RC207765], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PDE4A.




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-PDE4A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA501206], Dilution 1:50)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-PDE4A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA501206], Dilution 1:50)







Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-PDE4A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA501206], Dilution 1:50)



Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-PDE4A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA501206], Dilution 1:50)



Anti-PDE4A mouse monoclonal antibody ([TA501206]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PDE4A ([RC207765]).



HEK293T cells transfected with either [RC207765] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-PDE4A antibody ([TA501206]), and then analyzed by flow cytometry.