

Product datasheet for TA501270

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

IGF2BP2 Mouse Monoclonal Antibody [Clone ID: OTI1F9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1F9

Applications: FC, IF, IHC, WB

Recommended Dilution: WB 1:1000~2000, IHC 1:50, IF 1:100, FLOW 1:100

Reactivity: Human, Dog, Monkey, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human IGF2BP2(NP_006539) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 65.9 kDa

Gene Name: insulin like growth factor 2 mRNA binding protein 2

Database Link: NP 006539

Entrez Gene 303824 RatEntrez Gene 319765 MouseEntrez Gene 478662 DogEntrez Gene

700598 MonkeyEntrez Gene 10644 Human

Q9Y6M1

Background: This gene encodes a member of the IGF-II mRNA-binding protein (IMP) family. The protein

encoded by this gene contains several four KH domains and two RRM domains. It functions by binding to the 5' UTR of the insulin-like growth factor 2 (IGF2) mRNA and regulating IGF2 translation. Alternate transcriptional splice variants, encoding different isoforms, have been


characterized. [provided by RefSeq]






Synonyms: IMP-2; IMP2; VICKZ2


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY IGF2BP2 ([RC205673], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-IGF2BP2. Positive lysates [LY401961] (100ug) and [LC401961] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-IGF2BP2 monoclonal antibody.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-IGF2BP2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501270, Dilution 1:50)






Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-IGF2BP2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA501270, Dilution 1:50)




Anti-IGF2BP2 mouse monoclonal antibody (TA501270) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY IGF2BP2 ([RC205673]).



HEK293T cells transfected with either [RC205673] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-IGF2BP2 antibody (TA501270), and then analyzed by flow cytometry.





Flow cytometric Analysis of Jurkat cells, using anti-IGF2BP2 antibody (TA501270), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).