

## **Product datasheet for TA501614AM**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **NUMB Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1A7]**

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI1A7

Applications: FC, IF, WB

Recommended Dilution: WB 1:2000, IF 1:100, FLOW 1:100

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human NUMB (NP\_003735) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

Storage: Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 69.3 kDa

**Gene Name:** NUMB endocytic adaptor protein

Database Link: NP 003735

Entrez Gene 18222 MouseEntrez Gene 29419 RatEntrez Gene 8650 Human

P49757

**Background:** The protein encoded by this gene plays a role in the determination of cell fates during

development. The encoded protein, whose degradation is induced in a proteasome-dependent manner by MDM2, is a membrane-bound protein that has been shown to associate with EPS15, LNX1, and NOTCH1. Four transcript variants encoding different

isoforms have been found for this gene.








Synonyms: C14orf41; c14\_5527; S171

Protein Pathways: Notch signaling pathway

## **Product images:**








HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NUMB ([RC216230], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NUMB. Positive lysates [LY401229] (100ug) and [LC401229] (20ug) can be purchased separately from OriGene.

Western blot analysis of extracts (10ug) from 3 different cell lines by using anti-NUMB monoclonal antibody (1:200).

Anti-NUMB mouse monoclonal antibody ([TA501614]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY NUMB ([RC216230]).





HEK293T cells transfected with either [RC216230] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-NUMB antibody ([TA501614]), and then analyzed by flow cytometry.