

Product datasheet for TA501742BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

IVD Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1B10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1B10
Applications: IF, IHC, WB

Recommended Dilution: WB 1:500, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human, Dog, Rat, Monkey, Mouse

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human IVD (NP_002612) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 46.5 kDa

Gene Name: isovaleryl-CoA dehydrogenase

Database Link: NP 002216

Entrez Gene 24513 RatEntrez Gene 56357 MouseEntrez Gene 100856316 DogEntrez Gene

702867 MonkeyEntrez Gene 3712 Human

P26440

Background: Isovaleryl-CoA dehydrogenase (IVD) is a mitochondrial matrix enzyme that catalyzes the third

step in leucine catabolism. The genetic deficiency of IVD results in an accumulation of isovaleric acid, which is toxic to the central nervous system and leads to isovaleric acidemia. Alternatively spliced transcript variants encoding different isoforms have been found for this

gene. [provided by RefSeq]






Synonyms: ACAD2; IVDH


Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Valine, leucine and isoleucine degradation

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY IVD (Cat# [RC201077], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-IVD(Cat# [TA501742]). Positive lysates [LY432239] (100ug) and [LC432239] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-IVD monoclonal antibody.






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-IVD mouse monoclonal antibody ([TA501742]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-IVD mouse monoclonal antibody ([TA501742]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-IVD mouse monoclonal antibody ([TA501742]) at 1:150 dilution.






Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-IVD mouse monoclonal antibody ([TA501742]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-IVD mouse monoclonal antibody ([TA501742]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-IVD mouse monoclonal antibody ([TA501742]) at 1:150 dilution.






Immunohistochemical staining of paraffinembedded Human bladder tissue within the normal limits using anti-IVD mouse monoclonal antibody ([TA501742]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-IVD mouse monoclonal antibody ([TA501742]) at 1:150 dilution.



Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-IVD mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA501742], Dilution 1:50)





Anti-IVD mouse monoclonal antibody ([TA501742]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY IVD ([RC201077]).