

Product datasheet for TA501870S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ALDH1L1 Mouse Monoclonal Antibody [Clone ID: OTI5G8]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5G8

Applications: FC, IF, IHC, WB

Recommended Dilution: WB 1:500~2000, IF 1:100, FLOW 1:100

Reactivity: Human, Dog, Mouse, Rat

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ALDH1L1 (NP_036322) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 98.6 kDa

Gene Name: aldehyde dehydrogenase 1 family member L1

Database Link: NP 036322

Entrez Gene 64392 RatEntrez Gene 100855730 DogEntrez Gene 10840 Human

<u>075891</u>

Background: The protein encoded by this gene catalyzes the conversion of 10-formyltetrahydrofolate,

NADP, and water to tetrahydrofolate, NADPH, and carbon dioxide. The encoded protein belongs to the aldehyde dehydrogenase family and is responsible for formate oxidation in vivo. Deficiencies in this gene can result in an accumulation of formate and subsequent

methanol poisoning. [provided by RefSeq]






Synonyms: 10-fTHF; 10-FTHFDH; FDH; FTHFD


Protein Families: Druggable Genome

Protein Pathways: One carbon pool by folate


Product images:



Anti-ALDH1L1 mouse monoclonal antibody ([TA501870]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ALDH1L1 ([RC213720]).




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ALDH1L1 ([RC213720], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ALDH1L1. Positive lysates [LY415919] (100ug) and [LC415919] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-ALDH1L1 monoclonal antibody.






HEK293T cells transfected with either [RC213720] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ALDH1L1 antibody ([TA501870]), and then analyzed by flow cytometry.




Flow cytometric Analysis of Hela cells, using anti-ALDH1L1 antibody ([TA501870]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).




Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-ALDH1L1 mouse monoclonal antibody. ([TA501870])






Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-ALDH1L1 mouse monoclonal antibody. ([TA501870])




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-ALDH1L1 mouse monoclonal antibody. ([TA501870])




Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-ALDH1L1 mouse monoclonal antibody. ([TA501870])






Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-ALDH1L1 mouse monoclonal antibody. ([TA501870])



Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-ALDH1L1 mouse monoclonal antibody. ([TA501870])



Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-ALDH1L1 mouse monoclonal antibody. ([TA501870])