

## Product datasheet for TA501932M

### Glutathione Synthetase (GSS) Mouse Monoclonal Antibody [Clone ID: OT11A12]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OT11A12
Applications:	IHC, WB
Recommended Dilution:	WB: 1:500, IHC: 1:50-1:150
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human GSS (NP_000169) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.58 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	52.2 kDa
Gene Name:	glutathione synthetase
Database Link:	<a href="#">NP_000169</a> <a href="#">Entrez Gene 14854 Mouse</a> <a href="#">Entrez Gene 25458 Rat</a> <a href="#">Entrez Gene 2937 Human</a> <a href="#">P48637</a>
Background:	Glutathione is important for a variety of biological functions, including protection of cells from oxidative damage by free radicals, detoxification of xenobiotics, and membrane transport. The protein encoded by this gene functions as a homodimer to catalyze the second step of glutathione biosynthesis, which is the ATP-dependent conversion of gamma-L-glutamyl-L-cysteine to glutathione. Defects in this gene are a cause of glutathione synthetase deficiency. [provided by RefSeq]



[View online »](#)

**Synonyms:** GSHS; HEL-S-64p; HEL-S-88n


**Protein Families:** Druggable Genome

**Protein Pathways:** Glutathione metabolism, Metabolic pathways


### Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY GSS ([RC203174], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GSS. Positive lysates [LY424876] (100ug) and [LC424876] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffin-embedded Human colon tissue within the normal limits using anti-GSS mouse monoclonal antibody at 1:150 ([TA501932]). Heat-induced epitope retrieval by EDTA solution buffer pH 8.0 at 120°C for 3 min.



Immunohistochemical staining of paraffin-embedded Human liver tissue within the normal limits using anti-GSS mouse monoclonal antibody at 1:150 ([TA501932]). Heat-induced epitope retrieval by EDTA solution buffer pH 8.0 at 120°C for 3 min.