

Product datasheet for TA501953AM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SAMHD1 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1F6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1F6
Applications: FC, IF, WB

Recommended Dilution: WB 1:500~2000, IF 1:100, FLOW 1:100

Reactivity: Human, Monkey

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human SAMHD1 (NP_056289) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 72.0 kDa

Gene Name: SAM and HD domain containing deoxynucleoside triphosphate triphosphohydrolase 1

Database Link: NP 056289

Entrez Gene 709060 MonkeyEntrez Gene 25939 Human

Q9Y3Z3

Background: This gene may play a role in regulation of the innate immune response. The encoded protein

is upregulated in response to viral infection and may be involved in mediation of tumor necrosis factor-alpha proinflammatory responses. Mutations in this gene have been

associated with Aicardi-Goutieres syndrome. [provided by RefSeq]. COMPLETENESS: complete


on the 3' end.






Synonyms: CHBL2; DCIP; HDDC1; MOP-5; SBBI88


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SAMHD1 ([RC206013], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SAMHD1. Positive lysates [LY414526] (100ug) and [LC414526] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-SAMHD1 monoclonal antibody.




Western blot analysis of extracts (10ug) from a mouse cell line and 3 different mouse tissues by using anti-SAMHD1 monoclonal antibody (1:200).





Western blot analysis of extracts (10ug) from 10 Human tissue by using anti-SAMHD1 monoclonal antibody at 1:200 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: colon;10: spleen).



Anti-SAMHD1 mouse monoclonal antibody ([TA501953]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY SAMHD1 ([RC206013]).







Figure from citation: Immunofluorescence of SAMHD1 protein level by using anti-SAMHD1 antibody in human CD4+ T cells. <u>View Citation</u>







Flow cytometric Analysis of Hela cells, using anti-SAMHD1 antibody ([TA501953]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).