

Product datasheet for TA502030AM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

HDAC6 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI4C5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4C5

Applications: FC, IF, IHC, WB

Recommended Dilution: WB 1:1000~2000, IF 1:100, FLOW 1:100, IHC 1:150

Reactivity: Human, Dog, Rat, Monkey

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human HDAC6 (NP_006035) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 131.2 kDa

Gene Name: histone deacetylase 6

Database Link: NP 006035

Entrez Gene 84581 RatEntrez Gene 480907 DogEntrez Gene 711363 MonkeyEntrez Gene

<u>10013 Human</u>

Q9UBN7






Background: Histones play a critical role in transcriptional regulation, cell cycle progression, and

developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to class II of the histone deacetylase/acuc/apha family. It contains an internal duplication of two catalytic domains which appear to function independently of each other. This protein possesses histone deacetylase activity and represses transcription. [provided by RefSeq]


Synonyms: CPBHM; HD6; JM21; PPP1R90

Protein Families: Druggable Genome, Transcription Factors

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY HDAC6 (Cat# [RC209649], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HDAC6(Cat# [TA502030]). Positive lysates [LY401822] (100ug) and [LC401822] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-HDAC6 monoclonal antibody.






Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-HDAC6 mouse monoclonal antibody ([TA502030]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-HDAC6 mouse monoclonal antibody ([TA502030]) at 1:150 dilution.




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-HDAC6 mouse monoclonal antibody ([TA502030]) at 1:150 dilution.






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-HDAC6 mouse monoclonal antibody ([TA502030]) at 1:150 dilution.




Anti-HDAC6 mouse monoclonal antibody ([TA502030]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY HDAC6 ([RC209649]).




HEK293T cells transfected with either [RC209649] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-HDAC6 antibody ([TA502030]), and then analyzed by flow cytometry.





Flow cytometric Analysis of Hela cells, using anti-HDAC6 antibody ([TA502030]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-HDAC6 antibody ([TA502030]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).