

Product datasheet for **TA502197S**


RTL8B Mouse Monoclonal Antibody [Clone ID: OTI2D4]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI2D4
Applications:	IF, WB
Recommended Dilution:	WB 1:500, IF 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FAM127C (NP_001071641) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.4 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	13.1 kDa
Gene Name:	retrotransposon Gag like 8B
Database Link:	NP_001071641 Entrez Gene 441518 Human Q17RB0
Synonyms:	CXX1c; MAR8B



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FAM127C ([RC218852], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FAM127C. Positive lysates [LY421496] (100ug) and [LC421496] (20ug) can be purchased separately from OriGene.



Anti-FAM127C mouse monoclonal antibody ([TA502197]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FAM127C ([RC218852]).