

## Product datasheet for **TA502290M**


### ACY3 Mouse Monoclonal Antibody [Clone ID: OTI2G4]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI2G4
Applications:	FC, WB
Recommended Dilution:	WB 1:500, FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ACY3 (NP_542389) produced in 293T Cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.42 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	35.1 kDa
Gene Name:	aminoacylase 3
Database Link:	<a href="#">NP_542389</a> <a href="#">Entrez Gene 91703 Human</a> <a href="#">Q96HD9</a>
Synonyms:	ACY-3; HCBP1
Protein Pathways:	Alanine, aspartate and glutamate metabolism, Histidine metabolism

[View online »](#)


## Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ACY3 ([RC202287], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ACY3. Positive lysates [LY408962] (100ug) and [LC408962] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC202287] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ACY3 antibody ([TA502290]), and then analyzed by flow cytometry.



Flow cytometric Analysis of HeLa cells, using anti-ACY3 antibody ([TA502290]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-ACY3 antibody ([TA502290]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).