

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502388AM

ATP6V1F Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI5F10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5F10
Applications:	FC, IF, IHC, WB
Recommended Dilution:	WB 1:1000~2000, IHC 1:150, IF 1:100, FLOW 1:100
Reactivity:	Human, Dog, Rat, Monkey, Mouse
Host:	Mouse
lsotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human ATP6V1F (NP_004222) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	13.2 kDa
Gene Name:	ATPase H+ transporting V1 subunit F
Database Link:	<u>NP 004222</u> Entrez Gene 116664 RatEntrez Gene 475199 DogEntrez Gene 702112 MonkeyEntrez Gene 9296 Human Q16864




CRIGENE ATP6V1F Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI5F10] – TA502388AM


Background:This gene encodes a component of vacuolar ATPase (V-ATPase), a multisubunit enzyme that
mediates acidification of eukaryotic intracellular organelles. V-ATPase dependent organelle
acidification is necessary for such intracellular processes as protein sorting, zymogen
activation, receptor-mediated endocytosis, and synaptic vesicle proton gradient generation.
V-ATPase is composed of a cytosolic V1 domain and a transmembrane V0 domain. The V1
domain consists of three A and three B subunits, two G subunits plus the C, D, E, F, and H
subunits. The V1 domain contains the ATP catalytic site. The V0 domain consists of five
different subunits: a, c, c', c', and d. Additional isoforms of many of the V1 and V0 subunit
proteins are encoded by multiple genes or alternatively spliced transcript variants. This
encoded protein is the V1 domain F subunit protein. [provided by RefSeq, Jul 2008]

Synonyms:ATP6S14; VATF; Vma7Protein Pathways:Epithelial cell signaling in Helicobacter pylori infection, Metabolic pathways, Oxidative
phosphorylation, Vibrio cholerae infection


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ATP6V1F ([RC210728], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ATP6V1F. Positive lysates [LY418132] (100ug) and [LC418132] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-ATP6V1F monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).



Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-ATP6V1F mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA502388])

Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-ATP6V1F mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA502388])

Anti-ATP6V1F mouse monoclonal antibody ([TA502388]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ATP6V1F ([RC210728]).



HEK293T cells transfected with either [RC210728] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ATP6V1F antibody ([TA502388]), and then analyzed by flow cytometry.