

Product datasheet for TA502413BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Arginase 1 (ARG1) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4E6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4E6

Applications: FC, IHC, WB

Recommended Dilution: WB 1:1000~2000, IHC 1:50, FLOW 1:100

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ARG1 (NP_000036) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 34.6 kDa

Gene Name: arginase 1

Database Link: NP 000036

Entrez Gene 11846 MouseEntrez Gene 29221 RatEntrez Gene 383 Human

P05089

Background: Arginase catalyzes the hydrolysis of arginine to ornithine and urea. At least two isoforms of

mammalian arginase exist (types I and II) which differ in their tissue distribution, subcellular localization, immunologic crossreactivity and physiologic function. The type I isoform encoded by this gene, is a cytosolic enzyme and expressed predominantly in the liver as a component of the urea cycle. Inherited deficiency of this enzyme results in argininemia, an autosomal recessive disorder characterized by hyperammonemia. [provided by RefSeq]




Arginase 1 (ARG1) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4E6] – TA502413BM

Synonyms: arginase; arginase 1; liver; liver-type arginase; OTTHUMP00000017209; type I arginase


Protein Families: Druggable Genome

Protein Pathways: Arginine and proline metabolism, Metabolic pathways

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY ARG1 (Cat# [RC204649], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ARG1(Cat# [TA502413]). Positive lysates [LY424951] (100ug) and [LC424951] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-ARG1 monoclonal antibody.





Western blot analysis of extracts (10ug) from a mouse cell line and 3 different mouse tissues by using anti-ARG1 monoclonal antibody (1:200).






Figure A, Western blot analysis of overexpressed lysates(15ug per lane) from HEK293T cells transfected with empty plasmid ([PS100001], NC), human ARG2 plasmid ([RC206756], hARG2), human ARG1 plasmid ([RC204649], hARG1) using anti-ARG1 antibody [TA502413](1:500). Figure B, Western blot analysis of the same samples as figure A with anti-DDK antibody ([TA180144], 1:1000)



Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-ARG1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA502413])





HEK293T cells transfected with either [RC204649] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ARG1 antibody ([TA502413]), and then analyzed by flow cytometry.