

Product datasheet for **TA502537AM**


HAGHL Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI3H2]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI3H2
Applications:	FC, WB
Recommended Dilution:	WB 1:500, FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human HAGHL (NP_115680) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	31.1 kDa
Gene Name:	hydroxyacylglutathione hydrolase like
Database Link:	NP_115680 Entrez Gene 84264 Human Q6PII5
Synonyms:	MGC2605
Protein Pathways:	Pyruvate metabolism




[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY HAGHL (Cat# [RC200832], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HAGHL (Cat# [TA502537]). Positive lysates [LY410246] (100ug) and [LC410246] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC200832] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-HAGHL antibody ([TA502537]), and then analyzed by flow cytometry.



Flow cytometric Analysis of HeLa cells, using anti-HAGHL antibody ([TA502537]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-HAGHL antibody ([TA502537]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).