

Product datasheet for TA502568BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

NDUFB9 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI13H11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI13H11

Applications: FC, IHC, WB

Recommended Dilution: WB: 1:200 - 1:1000, IHC 1:150, FLOW 1:100

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 3-179 of human

NDUFB9(NP_0049906) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 21.7 kDa

Gene Name: NADH:ubiquinone oxidoreductase subunit B9

Database Link: NP 004996

Entrez Gene 66218 MouseEntrez Gene 299954 RatEntrez Gene 4715 Human

Q9Y6M9

Synonyms: B22; CI-B22; LYRM3; UQOR22


Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation,

Parkinson's disease







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY NDUFB9 (Cat# [RC200223], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NDUFB9 (Cat# [TA502568]). Positive lysates [LY417578] (100ug) and [LC417578] (20ug) can be purchased separately from OriGene.




Flow cytometric Analysis of Jurkat cells, using anti-NDUFB9 antibody ([TA502568]), (Red), compared to a nonspecific negative control antibody, (Blue).




Flow cytometric Analysis of Hela cells, using anti-NDUFB9 antibody ([TA502568]), (Red), compared to a nonspecific negative control antibody, (Blue).







Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-NDUFB9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA502568])




Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-NDUFB9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA502568])




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-NDUFB9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA502568])




Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-NDUFB9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA502568])




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-NDUFB9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA502568])




Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-NDUFB9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA502568])




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-NDUFB9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA502568])



Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-NDUFB9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA502568])



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-NDUFB9 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA502568])



Western blot analysis of extracts (10ug) from 6 different cell lines by using anti-NDUFB9 monoclonal antibody (1:200).