

## Product datasheet for **TA502595**


### **ANKRD53 Mouse Monoclonal Antibody [Clone ID: OTI2E4]**

#### **Product data:**


Product Type:	Primary Antibodies
Clone Name:	OTI2E4
Applications:	FC, WB
Recommended Dilution:	WB 1:500, FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-300 of human ANKRD53 (NP_079209) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.63 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	38 kDa
Gene Name:	ankyrin repeat domain 53
Database Link:	<a href="#">NP_079209</a> <a href="#">Entrez Gene 79998 Human</a> <a href="#">Q8N9V6</a>
Synonyms:	FLJ12056; FLJ36160




[View online »](#)

**Product images:**


HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ANKRD53 ([RC206701], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ANKRD53. Positive lysates [LY410969] (100ug) and [LC410969] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC206701] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-ANKRD53 antibody (TA502595), and then analyzed by flow cytometry.



Flow cytometric Analysis of HeLa cells, using anti-ANKRD53 antibody (TA502595), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-ANKRD53 antibody (TA502595), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).