

## Product datasheet for **TA502631S**


### UBXN2B Mouse Monoclonal Antibody [Clone ID: OTI1F7]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI1F7
Applications:	FC, IHC, WB
Recommended Dilution:	WB 1:2000, IHC 1:150, FLOW 1:100
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human UBXN2B (NP_001071087) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.41 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	36.9 kDa
Gene Name:	UBX domain protein 2B
Database Link:	<a href="#">NP_001071087</a> <a href="#">Entrez Gene 137886 Human</a> <a href="#">Q14CS0</a>
Synonyms:	p37




[View online »](#)

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY UBXN2B (Cat# [RC219137], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBXN2B(Cat# [TA502631]). Positive lysates [LY421462] (100ug) and [LC421462] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffin-embedded Human lymph node tissue within the normal limits using anti-UBXN2B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA502631])




Immunohistochemical staining of paraffin-embedded Human lymphoma tissue using anti-UBXN2B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA502631])



HEK293T cells transfected with either [RC219137] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-UBXN2B antibody ([TA502631]), and then analyzed by flow cytometry.



Flow cytometric Analysis of HeLa cells, using anti-UBXN2B antibody ([TA502631]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-UBXN2B antibody ([TA502631]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).