

Product datasheet for TA502643AM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Beclin 1 (BECN1) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1F1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1F1

Applications: FC, IF, IHC, WB

Recommended Dilution: WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human, Dog, Rat, Mouse

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human BECN1 (NP_003757) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Biotin

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 51.7 kDa **Gene Name:** beclin 1

Database Link: NP 003757

Entrez Gene 56208 MouseEntrez Gene 114558 RatEntrez Gene 480513 DogEntrez Gene 8678

<u>Human</u> Q14457

Background: Beclin-1 participates in the regulation of autophagy and has an important role in

development, tumorigenesis, and neurodegeneration (Zhong et al., 2009 [PubMed

19270693]). [supplied by OMIM]

Synonyms: ATG6; beclin1; VPS30






Beclin 1 (BECN1) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1F1] - TA502643AM


Protein Families: Druggable Genome

Protein Pathways: Regulation of autophagy


Product images:



Anti-BECN1 mouse monoclonal antibody ([TA502643]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY BECN1 ([RC201629]).




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-BECN1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA502643])




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY BECN1 (Cat# [RC201629], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BECN1(Cat# [TA502643]). Positive lysates [LY401239] (100ug) and [LC401239] (20ug) can be purchased separately from OriGene.






Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-BECN1 monoclonal antibody.





HEK293T cells transfected with either [RC201629] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-BECN1 antibody ([TA502643]), and then analyzed by flow cytometry.




Flow cytometric Analysis of Jurkat cells, using anti-BECN1 antibody ([TA502643]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).









Equivalent amounts of cell lysates (10 ug per lane) of wild-type HeLa cells (WT, Cat# LC810HELA) and BECN1-Knockout HeLa cells (KO, Cat# [LC834634]) were separated by SDS-PAGE and immunoblotted with anti-BECN1 monoclonal antibody [TA502643] (1:500). Then the blotted membrane was stripped and reprobed with anti-HSP90 antibody as a loading control.