

## Product datasheet for **TA502646S**

### UBXN2B Mouse Monoclonal Antibody [Clone ID: OTI1G8]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1G8
Applications:	FC, IF, WB
Recommended Dilution:	WB 1:2000, IF 1:100, FLOW 1:100
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human UBXN2B (NP_001071087) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.71 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	36.9 kDa
Gene Name:	UBX domain protein 2B
Database Link:	<a href="#">NP_001071087</a> <a href="#">Entrez Gene 137886 Human</a> <a href="#">Q14CS0</a>
Synonyms:	p37




[View online »](#)


## Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UBXN2B ([RC219137], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBXN2B. Positive lysates [LY421462] (100ug) and [LC421462] (20ug) can be purchased separately from OriGene.



Anti-UBXN2B mouse monoclonal antibody ([TA502646]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY UBXN2B ([RC219137]).



Flow cytometric Analysis of HeLa cells, using anti-UBXN2B antibody ([TA502646]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-UBXN2B antibody ([TA502646]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).