

## Product datasheet for TA502688BM

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### DCK Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3F5]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI3F5

**Applications:** FC, IF, IHC, WB

**Recommended Dilution:** WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DCK (NP\_000779) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 30.3 kDa

**Gene Name:** deoxycytidine kinase

Database Link: NP 000779

Entrez Gene 13178 MouseEntrez Gene 79127 RatEntrez Gene 1633 Human

P27707

**Background:** Deoxycytidine kinase (DCK) is required for the phosphorylation of several

deoxyribonucleosides and their nucleoside analogs. Deficiency of DCK is associated with resistance to antiviral and anticancer chemotherapeutic agents. Conversely, increased deoxycytidine kinase activity is associated with increased activation of these compounds to cytotoxic nucleoside triphosphate derivatives. DCK is clinically important because of its

relationship to drug resistance and sensitivity. [provided by RefSeq, Jul 2008]






**Synonyms:** MGC117410; MGC138632


**Protein Families:** Druggable Genome

**Protein Pathways:** Purine metabolism, Pyrimidine metabolism

# **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DCK ([RC210767], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DCK. Positive lysates [LY400272] (100ug) and [LC400272] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-DCK monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human) (1:200).







Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-DCK mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA502688])




Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-DCK mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA502688])




Anti-DCK mouse monoclonal antibody ([TA502688]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY DCK ([RC210767]).



HEK293T cells transfected with either [RC210767] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-DCK antibody ([TA502688]), and then analyzed by flow cytometry.



Flow cytometric Analysis of Hela cells, using anti-DCK antibody ([TA502688]), (Red), compared to a nonspecific negative control antibody, (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-DCK antibody ([TA502688]), (Red), compared to a nonspecific negative control antibody, (Blue).