

Product datasheet for **TA502699BM**

ECRG4 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3B10]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3B10
Applications:	FC, IF, WB
Recommended Dilution:	WB 1:500, IF 1:100, FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human C2orf40 (NP_115787) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	17 kDa
Gene Name:	ECRG4 augurin precursor
Database Link:	NP_115787 Entrez Gene 84417 Human Q9H1Z8
Synonyms:	ECRG4
Protein Families:	Secreted Protein, Transmembrane




[View online »](#)


Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY C2orf40 (Cat# [RC206239], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-C2orf40 (Cat# [TA502699]). Positive lysates [LY403168] (100ug) and [LC403168] (20ug) can be purchased separately from OriGene.



Flow cytometric Analysis of Jurkat cells, using anti-C2orf40 antibody ([TA502699]), (Red), compared to a nonspecific negative control antibody, (Blue).



Flow cytometric Analysis of HeLa cells, using anti-C2orf40 antibody ([TA502699]), (Red), compared to a nonspecific negative control antibody, (Blue).



Immunofluorescent staining of HeLa cells using anti-C2orf40 mouse monoclonal antibody ([TA502699]).