

## Product datasheet for **TA502707BM**


### **TBC1D4 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI5A6]**

#### **Product data:**


Product Type:	Primary Antibodies
Clone Name:	OTI5A6
Applications:	FC, IF, WB
Recommended Dilution:	WB 1:2000, IF 1:100, FLOW 1:100
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human TBC1D4 (NP_0055647) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	146.4 kDa
Gene Name:	TBC1 domain family member 4
Database Link:	<a href="#">NP_055647</a> <a href="#">Entrez Gene 210789 Mouse</a> <a href="#">Entrez Gene 306117 Rat</a> <a href="#">Entrez Gene 9882 Human</a> <a href="#">O60343</a>
Synonyms:	AS160; NIDDM5




[View online »](#)

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TBC1D4 ([RC212105], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TBC1D4. Positive lysates [LY415036] (100ug) and [LC415036] (20ug) can be purchased separately from OriGene.




Anti-TBC1D4 mouse monoclonal antibody ([TA502707]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY TBC1D4 ([RC212105]).



HEK293T cells transfected with either [RC212105] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-TBC1D4 antibody ([TA502707]), and then analyzed by flow cytometry.



Flow cytometric Analysis of Jurkat cells, using anti-TBC1D4 antibody ([TA502707]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).



Flow cytometric Analysis of HeLa cells, using anti-TBC1D4 antibody ([TA502707]), (Red), compared to a nonspecific negative control antibody (TA50011), (Blue).