

## Product datasheet for **TA502957**

### **NDUFAF7 Mouse Monoclonal Antibody [Clone ID: OTI6F10]**


#### **Product data:**

|                         |  |
|-------------------------|--|
| Product Type:           | Primary Antibodies   |
| Clone Name:             | OTI6F10  |
| Applications:           | FC, WB   |
| Recommended Dilution:   | WB 1:500, FLOW 1:100   |
| Reactivity:             | Human, Mouse   |
| Host:                   | Mouse  |
| Isotype:                | IgG2b  |
| Clonality:              | Monoclonal   |
| Immunogen:              | Full length human recombinant protein of human C2orf56 (NP_653337) produced in HEK293T cell.   |
| Formulation:            | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.   |
| Concentration:          | 1 mg/ml  |
| Purification:           | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)                              |
| Conjugation:            | Unconjugated   |
| Storage:                | Store at -20°C as received.  |
| Stability:              | Stable for 12 months from date of receipt.   |
| Predicted Protein Size: | 49.1 kDa   |
| Gene Name:              | NADH:ubiquinone oxidoreductase complex assembly factor 7   |
| Database Link:          | <a href="#">NP_653337</a><br><a href="#">Entrez Gene 73694 Mouse</a> <a href="#">Entrez Gene 55471 Human</a><br><a href="#">Q7L592</a> |
| Synonyms:               | C2orf56; MidA; PRO1853   |




[View online »](#)

## Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY C2orf56 ([RC202529], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-C2orf56. Positive lysates [LY403408] (100ug) and [LC403408] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC202529] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-C2orf56 antibody (TA502957), and then analyzed by flow cytometry.