

Product datasheet for **TA503858**

TRK fused gene (TFG) Mouse Monoclonal Antibody [Clone ID: OTI2C3]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI2C3
Applications:	IHC, WB
Recommended Dilution:	WB 1:2000, IHC 1:150
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human TFG(NP_001007566) produced in HEK293T cell.
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.45 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	43.3 kDa
Gene Name:	trafficking from ER to golgi regulator
Database Link:	NP_001007566 Entrez Gene 360709 Rat Entrez Gene 10342 Human Q92734
Background:	There are several documented fusion oncoproteins encoded partially by this gene. This gene also participates in several oncogenic rearrangements resulting in anaplastic lymphoma and mixoid chondrosarcoma, and may play a role in the NF-kappaB pathway. Multiple transcript variants have been found for this gene. [provided by RefSeq]
Synonyms:	HMSNP; SPG57; TF6; TRKT3




[View online »](#)

Protein Pathways: Pathways in cancer, Thyroid cancer


Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY TFG (Cat# [RC200293], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TFG (Cat# TA503858). Positive lysates [LY423514] (100ug) and [LC423514] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 1 cell line by using anti-TFG monoclonal antibody (1:200).




Equivalent amounts of cell lysates (10 ug per lane) of wild-type HeLa cells (WT, Cat# LC810HELA) and TFG-Knockout HeLa cells (KO, Cat# [LC831921]) were separated by SDS-PAGE and immunoblotted with anti-TFG monoclonal antibody TA503858 (1:200). Then the blotted membrane was stripped and reprobed with anti-PCNA antibody as a loading control.




Immunohistochemical staining of paraffin-embedded Human breast tissue within the normal limits using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)




Immunohistochemical staining of paraffin-embedded Human Kidney tissue within the normal limits using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)




Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human ovary tissue using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)




Immunohistochemical staining of paraffin-embedded Human pancreas tissue within the normal limits using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)




Immunohistochemical staining of paraffin-embedded Human endometrium tissue within the normal limits using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)




Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)




Immunohistochemical staining of paraffin-embedded Human prostate tissue within the normal limits using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)



Immunohistochemical staining of paraffin-embedded Carcinoma of Human bladder tissue using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)



Immunohistochemical staining of paraffin-embedded Human lymphoma tissue using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)



Immunohistochemical staining of paraffin-embedded Human tonsil within the normal limits using anti-TFG mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA503858)