

Product datasheet for **TA504256AM**

Frataxin (FXN) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1F2]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI1F2
Applications:	FC, WB
Recommended Dilution:	WB 1:2000, FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FXN(NP_000135) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	18.8 kDa
Gene Name:	frataxin
Database Link:	NP_000135 Entrez Gene 2395 Human Q16595
Background:	This nuclear gene encodes a mitochondrial protein which belongs to FRATAXIN family. The protein functions in regulating mitochondrial iron transport and respiration. The expansion of intronic trinucleotide repeat GAA results in Friedreich ataxia. Alternative splicing results in multiple transcript variants. [provided by RefSeq]
Synonyms:	CyaY; FA; FARR; FRDA; X25




[View online »](#)

Protein Families: Druggable Genome


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FXN ([RC204880], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FXN. Positive lysates [LY400054] (100ug) and [LC400054] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC204880] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-FXN antibody ([TA504256]), and then analyzed by flow cytometry.



Flow cytometric Analysis of Jurkat cells, using anti-FXN antibody ([TA504256]), (Red), compared to a nonspecific negative control antibody, (Blue).