

## Product datasheet for **TA504434M**

### **PEDF (SERPINF1) Mouse Monoclonal Antibody [Clone ID: OTI1F7]**

#### **Product data:**


Product Type:	Primary Antibodies
Clone Name:	OTI1F7
Applications:	WB
Recommended Dilution:	WB 1:1000
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human SERPINF1(NP_002606) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1.08 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	46.1 kDa
Gene Name:	serpin family F member 1
Database Link:	<a href="#">NP_002606</a> <a href="#">Entrez Gene 20317 Mouse</a> <a href="#">Entrez Gene 5176 Human</a> <a href="#">P36955</a>
Background:	The protein encoded by this gene is a member of the serpin family, although it does not display the serine protease inhibitory activity shown by many of the other serpin family members. The encoded protein is secreted and strongly inhibits angiogenesis. In addition, this protein is a neurotrophic factor involved in neuronal differentiation in retinoblastoma cells.


[View online »](#)

**Synonyms:** EPC-1; OI6; OI12; PEDF; PIG35

**Protein Families:** Druggable Genome, Secreted Protein

**Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY SERPINF1 (Cat# [RC204186], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SERPINF1(Cat# [TA504434]). Positive lysates [LY400927] (100ug) and [LC400927] (20ug) can be purchased separately from OriGene.