

## Product datasheet for **TA504547**


### **NANP Mouse Monoclonal Antibody [Clone ID: OTI4D5]**

#### **Product data:**


<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTI4D5
<b>Applications:</b>	IHC, WB
<b>Recommend Dilution:</b>	WB 1:500, IHC 1:150
<b>Reactivity:</b>	Human
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG2a
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length human recombinant protein of human NANP(NP_689880) produced in HEK293T cell.
<b>Formulation:</b>	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Concentration:</b>	1 mg/ml
<b>Purification:</b>	Purified from mouse ascites fluids by affinity chromatography
<b>Predicted Protein Size:</b>	27.6 kDa
<b>Gene Name:</b>	N-acetylneuraminic acid phosphatase
<b>Database Link:</b>	<a href="#">NP_689880</a> <a href="#">Entrez Gene 140838</a> <a href="#">Human</a>
<b>Synonyms:</b>	C20orf147; dj694B14.3; HDHD4
<b>Protein Pathways:</b>	Amino sugar and nucleotide sugar metabolism, Metabolic pathways



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NANP (RC205503, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NANP. Positive lysates LY407388 (100ug) and LC407388 (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffin-embedded Human lymph node tissue within the normal limits using anti-NANP mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA504547)