

Product datasheet for TA504565BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

HSD17B2 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4F10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4F10
Applications: FC, WB

Recommended Dilution: WB 1:2000, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human HSD17B2(NP_002144) produced in

HEK293T cell.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 42.6 kDa

Gene Name: hydroxysteroid 17-beta dehydrogenase 2

Database Link: NP 002144

Entrez Gene 3294 Human

P37059

Synonyms: EDH17B2; HSD17; SDR9C2


Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Androgen and estrogen metabolism, Metabolic pathways







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY HSD17B2 (Cat# [RC203293], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HSD17B2(Cat# [TA504565]). Positive lysates [LY419500] (100ug) and [LC419500] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC203293] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-HSD17B2 antibody ([TA504565]), and then analyzed by flow cytometry.



Flow cytometric Analysis of Jurkat cells, using anti-HSD17B2 antibody ([TA504565]), (Red), compared to a nonspecific negative control antibody, (Blue).