

Product datasheet for TA504604BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Ribophorin II (RPN2) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1G9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1G9
Applications: IHC, WB

Reactivity: WB 1:2000, IHC 1:150 **Reactivity:** Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human RPN2(NP_002942) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 67 kDa

Gene Name: ribophorin II

Database Link: NP 002942

Entrez Gene 6185 Human

P04844

Background: This gene encodes a type I integral membrane protein found only in the rough endoplasmic

reticulum. The encoded protein is part of an N-oligosaccharyl transferase complex that links high mannose oligosaccharides to asparagine residues found in the Asn-X-Ser/Thr consensus

motif of nascent polypeptide chains. This protein is similar in sequence to the yeast

oligosaccharyl transferase subunit SWP1. Alternatively spliced transcript variants encoding

different isoforms have been found for this gene. [provided by RefSeq, Oct 2008]






Ribophorin II (RPN2) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1G9] – TA504604BM

Synonyms: RIBIIR; RPN-II; RPNII; SWP1


Protein Families: Transmembrane

Protein Pathways: Metabolic pathways, N-Glycan biosynthesis


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY RPN2 (Cat# [RC205852], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RPN2(Cat# [TA504604]). Positive lysates [LY418998] (100ug) and [LC418998] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-RPN2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA504604])




Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-RPN2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA504604])





Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-RPN2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA504604])



Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-RPN2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA504604])