

## Product datasheet for **TA504616S**

### HSD17B2 Mouse Monoclonal Antibody [Clone ID: OTI3E9]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3E9
Applications:	FC, WB
Recommended Dilution:	WB 1:2000, FLOW 1:100
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human HSD17B2(NP_002144) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.88 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	42.6 kDa
Gene Name:	hydroxysteroid 17-beta dehydrogenase 2
Database Link:	<a href="#">NP_002144</a> <a href="#">Entrez Gene 3294 Human</a> <a href="#">P37059</a>
Synonyms:	EDH17B2; HSD17; SDR9C2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Androgen and estrogen metabolism, Metabolic pathways




[View online »](#)


Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HSD17B2 ([RC203293], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HSD17B2. Positive lysates [LY419500] (100ug) and [LC419500] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (10ug) from 1 cell line by using anti-HSD17B2 monoclonal antibody (1:200).



HEK293T cells transfected with either [RC203293] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-HSD17B2 antibody ([TA504616]), and then analyzed by flow cytometry.



Flow cytometric Analysis of Jurkat cells, using anti-HSD17B2 antibody ([TA504616]), (Red), compared to a nonspecific negative control antibody, (Blue).