

Product datasheet for TA504637S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

TDO2 Mouse Monoclonal Antibody [Clone ID: OTI1D1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1D1
Applications: WB

Recommended Dilution: WB 1:2000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human TDO2(NP_005642) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 47.7 kDa

Gene Name: tryptophan 2,3-dioxygenase

Database Link: NP 005642

Entrez Gene 56720 MouseEntrez Gene 64206 RatEntrez Gene 6999 Human

P48775

Background: Tryptophan 2,3-dioxygenase (EC 1.13.11.11) plays a role in catalyzing the first and rat-limiting

step in the kynurenine pathway, the major pathway of tryptophan metabolism. [supplied by

OMIM]


Synonyms: TDO; TO; TPH2; TRPO

Protein Pathways: Metabolic pathways, Tryptophan metabolism





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TDO2 ([RC202777], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TDO2. Positive lysates [LY401726] (100ug) and [LC401726] (20ug) can be purchased separately from OriGene.