

## Product datasheet for TA504879BM

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## TDO2 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI2A4]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI2A4

Applications: WB

Recommended Dilution: WB 1:2000

Reactivity: Human

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human TDO2(NP\_005642) produced in HEK293T

cell

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 47.7 kDa

**Gene Name:** tryptophan 2,3-dioxygenase

Database Link: NP 005642

Entrez Gene 6999 Human

P48775

**Background:** Tryptophan 2,3-dioxygenase (EC 1.13.11.11) plays a role in catalyzing the first and rate-

limiting step in the kynurenine pathway, the major pathway of tryptophan metabolism.

[supplied by OMIM]

**Synonyms:** TDO; TO; TPH2; TRPO

**Protein Pathways:** Metabolic pathways, Tryptophan metabolism





## **Product images:**




Figure A, Western blot analysis of overexpressed lysates(15ug per lane) from HEK293T cells transfected with empty plasmid ([PS100001], lane 1), human TDO2 plasmid ([RC206592], lane 2), mouse TDO2 plasmid ([MR206377], lane 3) using anti-TDO2 antibody [TA504879] (1:500). Figure B, Western blot analysis of the same samples as figure A with anti-DDK antibody ([TA180144], 1:1000)