

## **Product datasheet for TA505283BM**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### B3GNT1 (B3GNT2) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3A9]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI3A9

**Applications:** IF, IHC, WB

Recommended Dilution: WB 1:2000, IHC 1:150, IF 1:100

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human B3GNT2(NP\_006568) produced in HEK293T

cell.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 45.8 kDa

**Gene Name:** UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase 2

Database Link: NP 006568

Entrez Gene 53625 MouseEntrez Gene 305571 RatEntrez Gene 10678 Human

**09NY97** 

**Background:** This gene encodes a member of the beta-1,3-N-acetylglucosaminyltransferase family. This

enzyme is a type II transmembrane protein. It prefers the substrate of lacto-N-neotetraose, and is involved in the biosynthesis of poly-N-acetyllactosamine chains. [provided by RefSeq,

Jul 20081





# B3GNT1 (B3GNT2) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3A9] – TA505283BM

Synonyms: 3-Gn-T1; 3-Gn-T2; B3GN-T2; B3GNT-2; B3GNT-1; beta-1; beta-1;


BETA3GNT; BGnT-2

**Protein Families:** Transmembrane


**Protein Pathways:** Glycosphingolipid biosynthesis - lacto and neolacto series, Keratan sulfate biosynthesis,

Metabolic pathways

### **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY B3GNT2 (Cat# [RC208100], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-B3GNT2(Cat# [TA505283]). Positive lysates [LY416552] (100ug) and [LC416552] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-B3GNT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA505283])






Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-B3GNT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA505283])



Anti-B3GNT2 mouse monoclonal antibody ([TA505283]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY B3GNT2 ([RC208100]).



Immunofluorescent staining of HeLa cells using anti-B3GNT2 mouse monoclonal antibody ([TA505283]).