

Product datasheet for **TA505298AM**


NKIRAS1 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1D2]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI1D2
Applications:	IF, WB
Recommended Dilution:	WB 1:2000, IF 1:100
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human NKIRAS1(NP_065078) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	21.5 kDa
Gene Name:	NFKB inhibitor interacting Ras like 1
Database Link:	NP_065078 Entrez Gene 69721 Mouse Entrez Gene 305751 Rat Entrez Gene 28512 Human Q9NYS0
Synonyms:	kappaB-Ras1; KBRAS1



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NKIRAS1 ([RC204248], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NKIRAS1. Positive lysates [LY412541] (100ug) and [LC412541] (20ug) can be purchased separately from OriGene.



Anti-NKIRAS1 mouse monoclonal antibody ([TA505298]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY NKIRAS1 ([RC204248]).