

# **Product datasheet for TA505353S**

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **DEF6 Mouse Monoclonal Antibody [Clone ID: OTI4B5]**

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4B5
Applications: IHC, WB

**Reactivity:** WB 1:2000, IHC 1:150 **Reactivity:** Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DEF6(NP\_071330) produced in HEK293T

cell.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Predicted Protein Size:** 73.7 kDa

**Gene Name:** DEF6 guanine nucleotide exchange factor

Database Link: NP 071330

Entrez Gene 23853 MouseEntrez Gene 309642 RatEntrez Gene 50619 Human

Q9H4E7

Background: DEF6, or IBP, is a guanine nucleotide exchange factor (GEF) for RAC (MIM 602048) and CDC42

(MIM 116952) that is highly expressed in B and T cells (Gupta et al., 2003 [PubMed

12923183]). [supplied by OMIM]


Synonyms: IBP; SLAT; SWAP70L


Protein Families: Druggable Genome






# **Product images:**







Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-DEF6 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, [TA505353])