

Product datasheet for **TA505462BM**


FNDC4 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI3B11]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI3B11
Applications:	IF, WB
Recommended Dilution:	WB 1:2000, IF 1:100
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FNDC4(NP_073734) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	25 kDa
Gene Name:	fibronectin type III domain containing 4
Database Link:	NP_073734 Entrez Gene 64339 Mouse Entrez Gene 64838 Human Q9H6D8
Synonyms:	FRCP1
Protein Families:	Transmembrane



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY FNDC4 (Cat# [RC207471], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FNDC4 (Cat# [TA505462]). Positive lysates [LY411537] (100ug) and [LC411537] (20ug) can be purchased separately from OriGene.



Anti-FNDC4 mouse monoclonal antibody ([TA505462]) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FNDC4 ([RC207471]).