

## Product datasheet for **TA505698AM**


### TRMT1 (TRMU) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI1C6]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1C6
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human TRMU (NP_060476) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	47.6 kDa
Gene Name:	tRNA mitochondrial 2-thiouridylase
Database Link:	<a href="#">NP_060476</a> <a href="#">Entrez Gene 72026 Mouse</a> <a href="#">Entrez Gene 362976 Rat</a> <a href="#">Entrez Gene 55687 Human</a> <a href="#">O75648</a>
Background:	This gene is a member of the trmU family. It encodes a mitochondria-specific tRNA-modifying enzyme that is required for the 2-thio modification of 5-taurinomethyl-2-thiouridine tRNA-Lys on the wobble position of the anticodon. [provided by RefSeq, Jul]
Synonyms:	LCAL3; MTO2; MTU1; TRMT; TRMT1; TRNT1



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TRMU ([RC219032], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRMU. Positive lysates [LY402642] (100ug) and [LC402642] (20ug) can be purchased separately from OriGene.