

Product datasheet for TA505835BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

NRAS Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI5G7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5G7
Applications: IHC, WB

Reactivity: WB 1:2000, IHC 1:150 **Reactivity:** Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NRAS (NP_002515) produced in HEK293T

cell

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 21 kDa

Gene Name: NRAS proto-oncogene, GTPase

Database Link: NP 002515

Entrez Gene 18176 MouseEntrez Gene 24605 RatEntrez Gene 4893 Human

P01111





Background:

This is an N-ras oncogene encoding a membrane protein that shuttles between the Golgi apparatus and the plasma membrane. This shuttling is regulated through palmitoylation and depalmitoylation by the ZDHHC9-GOLGA7 complex. The encoded protein, which has intrinsic GTPase activity, is activated by a guanine nucleotide-exchange factor and inactivated by a GTPase activating protein. Mutations in this gene have been associated with somatic rectal cancer, follicular thyroid cancer, autoimmune lymphoproliferative syndrome, Noonan syndrome, and juvenile myelomonocytic leukemia. [provided by RefSeq, Jun 2011]


Synonyms: ALPS4; CMNS; N-ras; NCMS; NRAS1; NS6

Protein Families: Druggable Genome

Protein Pathways: Acute myeloid leuk


Acute myeloid leukemia, Axon guidance, B cell receptor signaling pathway, Bladder cancer, Chemokine signaling pathway, Chronic myeloid leukemia, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Gap junction, Glioma, GnRH signaling pathway, Insulin signaling pathway, Long-term depression, Long-term potentiation, MAPK signaling pathway, Melanogenesis, Melanoma, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Non-small cell lung cancer, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton, Renal cell carcinoma, T cell receptor signaling pathway, Thyroid cancer, Tight junction, VEGF signaling pathway

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY NRAS (Cat# [RC202681], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NRAS(Cat# [TA505835]). Positive lysates [LY400901] (100ug) and [LC400901] (20ug) can be purchased separately from OriGene.






Western blot analysis of anti-NRAS monoclonal antibodiest, TA505835. Incubation: 1:500, 1h.

1: lysate of 293T transfected with HRAS plasmid, RC225202

2: lysate of 293T transfected with NRAS plasmid, RC202681

3. lysate of 293T transfected with KRAS plasmid, RC222697



HEK293T cells were transfected with the 3 different overexpression plasmids (1:HRAS, Cat# [RC225202];2: NRAS, Cat# [RC202681]; 3:KRAS, Cat# [RC222697]) for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-flag antibody (Cat# [TA180144], 1:1000) or anti-NRAS mouse monoclonal antibody (Cat# [TA505835], 1:500).

Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-NRAS mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA505835])